

Travel PAP Comparison Chart

React Health Luna TravelPAP



ResMed AirMini™



Philips DreamstationGo™ (Currently Recalled and unavailable for purchase)



Transcend Micro



| | | | | |
|--------------------------|------------------------------|--|---|---|
| Weight | .88 lbs (14.1 oz) | .66 lbs (10.6 oz) | 1.88 lbs (29.8 oz) | .48 lbs (7.7 oz) |
| Dimensions | 6.2"L X 2.6"W x 2.8"H | 5.4" W x 3.3" D x 2.1" H | 5.9" L x 5.9" W x 2.3" H | 3.6"L x 3.6"W x 2.4"H |
| Modes | CPAP, Auto | CPAP, AutoSet (with Response option) and AutoSet for Her | CPAP, APAP | Auto |
| Filter | Disposable | Disposable | Reusable/ Disposable | Reusable |
| Operating pressure range | 4 to 20 cm H2O | 4 to 20 cm H2O | 4 to 20 cm H2O | 4 to 20 cm H2O |
| Noise Level | < 30 dBA | 29 dBA | 30 dBA | 27dBa |
| Humidification | N/A | HumidX, HumidX Plus, and HumidX F20 | Optional Heated Humidifier | Airmist™ Heat and Moisture exchange (HME) Adapter |
| Elevation Operation | Sea level to 7,500 feet | Sea level to 8,500 feet | Sea Level to 7,500 feet | Sea Level to 8,000 feet |
| Connectivity/App | LightTrip (US) APP | AirMini APP | DreamMapper | MySleepDash APP |
| Ramp Time | 0 to 60 mins | 5 to 45 min | 0 to 45 min | 1-to 45 mins |
| Masks Required? | No, Standard Mask | Yes | No, Standard Mask | No, Standard Mask |
| Can be used on airplane | Yes (must have power source) | Yes (must have power source) | Yes, power source or optional external battery pack | Yes, power source or optional external battery pack |
| Warranty | 2 Years | 2 years | 2 Years | 2 Years |
| External Battery | No | No | Optional Battery Pack | Optional Battery Pack |
| Bluetooth Connectivity | Yes | Yes | Yes | Yes |
| DC Adapter / Charger | Yes | Yes | No | No |

*Micro, EZEX and Airmist are a trademarks of Transcend.

*AirMini is a trademark of ResMed, Inc.

*Dreamstation and DreamMapper are trademarks of Koninklijke Philips N.V.