



React Health's best selling full face mask offers the following features for patient comfort and compliance:

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for excellent seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- Swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet nights sleep



REACTHEALTH

Siesta
Full Face Mask



The Siesta Full Face Mask is labeled on the back of the mask.

Allows for ease of identification for troubleshooting and resupply



Mask Part Numbers

- SFF1000 - Siesta Full Face Mask - All Sizes Kit
- SFF1001 - Siesta Full Face Mask - Small
- SFF1002 - Siesta Full Face Mask - Medium
- SFF1003 - Siesta Full Face Mask - Large

Replacement Supplies Part Numbers

- SFF3001 - Replacement Cushion - Small
- SFF3002 - Replacement Cushion - Medium
- SFF3003 - Replacement Cushion - Large
- SFF2001 - Replacement Headgear
- SFF4001- Frame w cushion/pillow- Small
- SFF4002- Frame w cushion/pillow- Medium
- SFF4003- Frame w cushion/pillow- Large



REACTHEALTH

www.reacthealth.com
info@reacthealth.com
Phone: 863-226-6285