

# VC Pro

## NON-INVASIVE VENTILATOR



## WITH COUGH-ASSIST

- ▶ Provides invasive and non-invasive ventilation, integrated high flow therapy, and cough-assist
- ▶ Delivers a comprehensive set of modes and settings including pressure and volume control, as well as advanced leak compensation
- ▶ Touch Button Cough™ is activated with the touch of a button to make airway clearance easy
- ▶ One patient circuit provides both ventilation and cough-assist, supporting patient breathing throughout cough therapy

**Contact your React Health/Ventec representative:  
ventec-info@reacthealth.com or call 844-698-6276**



Invasive, and  
Non-invasive

# HFT

High Flow Therapy



Cough Assist



Up to 9 hours of  
battery life



Hospital to Home




Pediatric to Adult



## VC+Pro Offers 9 Ventilation Modes + Cough-Assist

- Spontaneous
- Bi-Level (with functionality similar to S/T, Timed, and BiLevel ventilation)
- Assist/Control-Pressure
- Assist/Control-Volume
- SIMV-Pressure (including CPAP)
- SIMV-Volume
- Vol. Targeted-PS (with functionality similar to AVAPS® and iVAPS®)
- Vol. Targeted-PC (with functionality similar to PRVC)
- Vol. Targeted-SIMV (with functionality similar to SIMV+PRVC)

### Patient Breathing Supported Throughout Cough Therapy



Configuration	Ventilation	Cough	High Flow Nasal Cannula Therapy	High Pressure External O2 and FiO2 Monitor	Low Pressure External Oxygen
<b>V+C+Pro</b> (PRT-00973-000)	✓	✓	✓	✓	✓
<b>V+Pro</b> (PRT-01198-000)	✓	—	✓	✓	✓