



REACT HEALTH

## New Luna G3 CPAP and Auto-CPAP



Luna G3 CPAP  
SKU: LG3500

Luna G3 Auto-CPAP  
SKU: LG3600

React Health's all new Luna G3 CPAP and Auto-CPAP devices come standard with all the features that make sense for you and your patients. Our standard integrated cellular modem makes capturing patient compliance a breeze.

Our standard patient-forward features include integrated heated tubing, pre-heat humidifier function and Auto Ramp offering added comfort to patient therapy as well as a bedside clock and resupply reminders for the convenience of the patient. You will find these features and more, in a smaller footprint, with an easy-to-read full color screen and intuitive interface.

# Device Specifications

Device Dimensions (L x W x H)	10.4" x 5.7" x 4.5"
Weight	3.75 lbs.
Mode of Operations	CPAP (LG3500), CPAP, Auto-CPAP (LG3600)
SD Card	Standard
Pressure Range	4 to 20cmH2O (in 0.5 cm increments)
Ramp Time	0 to 60 minutes
Sound	<26 dBA when device is at 10cmH2O
Water Capacity	350 ml (to maximum fill line)

## Replacement Parts

LG3500	Luna G3 CPAP Device - HT
LG3600	Luna G3 APAP Device - HT
LG34510	Luna G3 Replacement Water Chamber
LG3HT	Luna G3 Integrated Heated Tubing
UN100C	3B Ultra Noir CPAP Standard Tubing (Non-Heated) cs of 50
LG31020BX*	Luna G3 Disposable Air Filter - 6 Filters
LG38008	Luna G3 Reusable Filter - Box of 12 Filters (Black)
LG3B	Luna G3 Carry Bag
LPG3	Luna G3 Power Supply Adapter
LPAC	Luna G2/G3 Power Cord
PRTG301	Luna G3 Repl Filter Cover- Black
PRTG302	Luna G3 Repl Filter Cover- Ivory



React Health offers more ways to connect



*\*Additional Fees may apply for Wi-Fi and Cellular*



**REACTHEALTH**

[www.reacthealth.com](http://www.reacthealth.com)

[info@reacthealth.com](mailto:info@reacthealth.com)

\* Luna G3 has interchangeable filters. You may choose either the LG-31020BX (A7038) or LG38008BX (A7039) as a filter for the device.