

VC Pro

NON-INVASIVE VENTILATOR



WITH COUGH-ASSIST

- ▶ Provides invasive and non-invasive ventilation, integrated high flow therapy, and cough-assist
- ▶ Delivers a comprehensive set of modes and settings including pressure and volume control, as well as advanced leak compensation
- ▶ Touch Button Cough™ is activated with the touch of a button to make airway clearance easy
- ▶ One patient circuit provides both ventilation and cough-assist, supporting patient breathing throughout cough therapy

**Contact your React Health/Ventec representative:
ventec-info@reacthealth.com or call 844-698-6276**



Invasive, and
Non-invasive

HFT

High Flow Therapy



Cough Assist



Up to 9 hours of
battery life



Hospital to Home



Pediatric to Adult




VC+Pro Offers 9 Ventilation Modes + Cough-Assist

- Spontaneous
- BiLevel (with functionality similar to S/T, Timed, and BiLevel ventilation)
- Assist/Control-Pressure
- Assist/Control-Volume
- SIMV-Pressure (including CPAP)
- SIMV-Volume
- Vol. Targeted-PS (with functionality similar to AVAPS® and iVAPS®*)
- Vol. Targeted-PC (with functionality similar to PRVC)
- Vol. Targeted-SIMV (with functionality similar to SIMV+PRVC)

* - AVAPS® and iVAPS® are registered trademarks of Philips RS North America LLC

Patient Breathing Supported Throughout Cough Therapy



Configuration	Ventilation	Cough	High Flow Nasal Cannula Therapy	High Pressure External O2 and FiO2 Monitor	Low Pressure External Oxygen
V+C+Pro (PRT-00973-000)	✓	✓	✓	✓	✓
V+Pro (PRT-01198-000)	✓	—	✓	✓	✓