



REACTHEALTH

BILEVEL COMPARISON CHART

MANUFACTURER

React Health

Philips

ResMed



BRAND	LUNA G3 BILEVEL 25A	DREAMSTATION® BIPAP PRO*	AIRCURVE® 10 VAUTO
Dimensions	5.70 x 10.48 x 4.49 in.	11.7 x 7.6 x 3.3 in.	4.57 X 10.04 X 5.91 in.
Weight	3.74 lbs.	4.37 lbs.	2.75 lbs.
Humidifications	Standard	Fixed, Adaptive, Heated Tube modes	Standard, Integrated, Fixed or Adaptive
Heated Tubing	Standard	Optional	Optional
Connectivity	Code connect, QR Code, SD Card, Integrated Cellular Modem	Integrated Bluetooth, optional Wi-Fi & cellular modems	Wi-Fi, SD Card
SD Data Capacity	Minimum 12 months	1 year on-board: <6 months	365 sessions
Filters	Reusable inlet filter	Reusable pollen filter, disposable ultra-fine filter	Disposable Antibacterial Filter
Modes	CPAP, AutoCPAP, S, AutoS	Fixed CPAP, Fixed BiPAP	CPAP, BiLevel, Spontaneous
Rise Time	150-400 ms	100-400 ms	150-900 ms
Ramp Time	0 to 60 minutes	5 to 45 minutes (5 minute increments)	5 to 45 minutes, user adjustable
Ventilation Pressure	4 to 25 cm H2O	4 to 25 cm H2O	3-25cm H2O
Weight Restrictions	Over 66 lbs	Over 66 lbs.	Over 66 lbs.

*The Dreamstation BiPAP is not currently available for sale
 Dreamstation® is a registered trademark of Koninklijke Philips N.V. (Royal Philips)
 Philips AirCurve® is a registered trademark of ResMed family of companies.